


U 504
Basket
For the optimum loading of cutlery.



- With 9 plastic cutlery holders
- Compartment size 106 x 106 mm

EAN: 4002515250112 / Material number: 09508670 / Old Material Number: 67650401D

Product category	
Fresh-water dishwashers	•
Tank dishwashers	•
Throughfeed tank dishwashers	•
Machine affiliation	
G 8066	•
PG 8041	•
PG 8043	•
PG 8045	•
PG 8055	•
PG 8055 U	•
PG 8056	•
PG 8056 U	•
PG 8057 TD	•
PG 8057 TD U	•
PG 8058	•
PG 8058 U	•
PG 8059	•
PG 8059 U	•
PG 8060	•
PG 8061	•
PG 8061 U	•
PG 8165	•
PG 8166	•
PG 8169	•
PG 8172	•
PG 8172 Eco	•
PTD 901	•
PG 8096	•
PG 8096 U	•
PG 8099	•
PG 8099 U	•
PFD 400	•
PFD 400 U	•
PFD 401	•
PFD 401 U	•
PFD 402	•
PFD 402 U	•
PFD 404	•
PFD 404 U	•
PFD 405	•
PFD 405 U	•
PFD 407	•
PFD 407 U	•
PFD 408 U	•
PTD 702	•
PTD 703	•
PTD 704	•
Product category	
Basket for cutlery	•

U 504

Basket

For the optimum loading of cutlery.



EAN: 4002515250112 / Material number: 09508670 / Old Material Number: 67650401D

Product characteristics	
Material	Wire, rilsan-coated, Plastic
Color	Grey, White
Number of compartments	9
Compartment width in inches	4 3/16 (106)
Compartment length in inches	4 3/16 (106)
Dimensions and weight	
External dimensions, net height in inches	5 3/4 (145)
External dimensions, net width in inches	19 11/16 (500)
External dimensions, net depth in inches	19 11/16 (500)
External dimensions, gross height in inches	10 1/4 (260)
External dimensions, gross width in inches	20 1/2 (520)
External dimensions, gross depth in inches	20 7/8 (530)
Net weight in lbs	5 (2.24)
Gross weight in lbs	8 (3.21)