


E 167
Cutlery basket
For holding cutlery.



- Suitable for plate basket U 506
- 7 Compartments

EAN: 4002515250167 / Material number: 09508730 / Old Material Number: 67616701D

Product category	
Fresh-water dishwashers	•
Tank dishwashers	•
Throughfeed tank dishwashers	•
Machine affiliation	
G 8066	•
PG 8041	•
PG 8043	•
PG 8045	•
PG 8055	•
PG 8055 U	•
PG 8056	•
PG 8056 U	•
PG 8057 TD	•
PG 8057 TD U	•
PG 8058	•
PG 8058 U	•
PG 8059	•
PG 8059 U	•
PG 8060	•
PG 8061	•
PG 8061 U	•
PG 8165	•
PG 8166	•
PG 8169	•
PG 8172	•
PG 8172 Eco	•
PTD 901	•
PG 8096	•
PG 8096 U	•
PG 8099	•
PG 8099 U	•
PFD 400	•
PFD 400 U	•
PFD 401	•
PFD 401 U	•
PFD 402	•
PFD 402 U	•
PFD 404	•
PFD 404 U	•
PFD 405	•
PFD 405 U	•
PFD 407	•
PFD 407 U	•
PFD 408 U	•
PTD 702	•
PTD 703	•
PTD 704	•
Product category	
Insert for cutlery	•

E 167
Cutlery basket
For holding cutlery.



EAN: 4002515250167 / Material number: 09508730 / Old Material Number: 67616701D

Product characteristics	
Material	Plastic
Color	White
Number of compartments	7
Dimensions and weight	
External dimensions, net height in inches	4 3/16 (106)
External dimensions, net width in inches	17 3/8 (440)
External dimensions, net depth in inches	2 (50)
External dimensions, gross height in inches	3 15/16 (100)
External dimensions, gross width in inches	5 15/16 (150)
External dimensions, gross depth in inches	29 9/16 (750)
Net weight in lbs	5 (2.24)
Gross weight in lbs	8 (3.21)