

HUBB 71

Genuine Miele multi-purpose tray with PerfectClean finish.

EAN: 4002515261286 / Material number: 09519840 / Old Material Number: 22996250D



- As a drip tray for meat juices/fat from roasting and grilling
- For use when roasting with the grilling and roasting insert
- Fits into FlexiClips or directly into side racks
- For H 2000/6000/7000 and DGC 6000/7000 appliances
- Dimensions: 18" x 15" x 1.4" (44.8 cm x 38.6 cm x 3.475 cm)
- Especially easy cleaning thanks to Perfect Clean finish