

E 900-6
Mobile unit
for the optimum cleaning of instruments on 6 levels.



- Equipped with connections for hollow instruments
- Loading capacity: 18 DIN mesh trays*
- Equipped with 3 integrated spray arms
- The loading areas for the individual levels are different
- Magnetic strip for mobile unit recognition (incl. magnets ML)

* DIN mesh tray 485 x 254 x 50 mm

EAN: 4002515263358 / Material number: 09525120 / Old Material Number: 69590007D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
PG 8527	•
PG 8528	•
Product category	
Mobile unit for inserts	•
Mobile unit for mesh trays	•
Mobile unit for MIS instruments	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Reprocessing of surgical instruments	•
Capacity	
DIN mesh trays [number]	18
SPRI mesh trays [number]	12
Lumened instruments [number]	16
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in inches	21 13/16 (553)
External dimensions, net width in inches	25 1/16 (636)
External dimensions, net depth in inches	31 1/8 (790)
External dimensions, gross height in inches	26 13/16 (680)
External dimensions, gross width in inches	22 1/16 (560)
External dimensions, gross depth in inches	25 5/8 (650)
Net weight in lbs	29 (12.85)
Gross weight in lbs	29 (12.85)
Accessory pack	
Magnetic strip for automatic mobile unit recognition (incl. magnets) [number]	1