

## Stand 1-80

### Conversion kit

For converting a built-under appliance to a free-standing appliance.



- Suitable for PG 8056 and PG 8061 dishwashers
- Contains toe kick with two brackets, lid and rear panel
- Available in stainless steel

EAN: 4002515426982 / Material number: 09800680 / Old Material Number: 67211603USA

#### Machine affiliation

PG 8055 U

•

PG 8056 U

•

PG 8057 TD U

•

PG 8058 U

•

PG 8059 U

•

PG 8061 U

•

#### Dimensions and weight

External dimensions, gross height in inches

7 7/8 (200)

External dimensions, gross width in inches

24 13/16 (630)

External dimensions, gross depth in inches

24 7/16 (620)

Net weight in lbs

29 (13.00)

Gross weight in lbs

31 (14.00)