


A 102

Upper basket

Height adjustable, for the optimum loading of inserts.



- Vertical clearance 205 +/- 30 mm, width 475 mm, depth 443 mm
- With integrated spray arm

EAN: 4002515460634 / Material number: 09862360 / Old Material Number: 69610200D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Thermal disinfectors	•
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8591	•
PG 8592	•
PG 8593	•
<b>Product category</b>	
Carrier/basket carrier	•
Upper basket	•
<b>Product characteristics</b>	
Material	Stainless steel
Color	Stainless Steel
<b>Application area</b>	
Reprocessing of laboratory glassware	•
Reprocessing of GYN instruments	•
Reprocessing of ENT instruments	•
Reprocessing of dental instruments	•
Reprocessing of medical shoes	•
Reprocessing of basins	•
Reprocessing of kidney bowls	•
Reprocessing of surgical instruments	•
<b>Capacity</b>	
DIN mesh trays [number]	2
SPRI mesh trays [number]	1
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, net height in inches	8 1/8 (206)
External dimensions, net width in inches	20 13/16 (528)
External dimensions, net depth in inches	20 3/4 (527)
External dimensions, gross height in inches	9 1/16 (230)
External dimensions, gross width in inches	21 5/16 (540)
External dimensions, gross depth in inches	22 1/4 (565)
Net weight in lbs	7 (3.02)
Gross weight in lbs	9 (3.78)