

A 150

Lower basket

For the optimum loading of injector modules or inserts.



- Lower basket with two module connections
- for holding up to 2 injector modules or inserts
- Automatic closing of the connections when not in use
- Clearances: width 19 4/9" (494 mm), length 17 5/16" (440 mm)

EAN: 4002515460672 / Material number: 09862400 / Old Material Number: 69615000D

Product category	
Washer-disinfectors, laboratory	•
Machine affiliation	
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
Product category	
Basket with module connection	•
Lower basket	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black
Application area	
Reprocessing of pipettes/vials	•
Reprocessing of laboratory glassware	•
Reprocessing of butyrometers	•
Reprocessing of measuring cylinders	•
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in inches	6 1/8 (154)
External dimensions, net width in inches	20 7/8 (529)
External dimensions, net depth in inches	21 1/2 (546)
External dimensions, gross height in inches	6 15/16 (175)
External dimensions, gross width in inches	21 1/2 (545)
External dimensions, gross depth in inches	22 1/4 (565)
Net weight in lbs	6 (2.61)
Gross weight in lbs	8 (3.29)