


A 200  
Mobile unit  
for the optimum loading of 38 pipettes.



- for 10 100 ml pipettes in row 1, distance 20 mm
- for 14 25 ml pipettes in a 2nd row, distance 26 mm
- for 14 10 ml pipettes in a 3rd row, distance 26 mm

EAN: 4002515460696 / Material number: 09862420 / Old Material Number: 69620000D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PG 8583	•
PG 8583 CD	•
PG 8593	•
<b>Product category</b>	
Mobile unit for pipettes	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black
<b>Application area</b>	
Reprocessing of pipettes/vials	•
<b>Capacity</b>	
Pipettes [number]	38
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in inches	15 11/16 (397)
External dimensions, net width in inches	20 7/8 (529)
External dimensions, net depth in inches	21 1/2 (546)
External dimensions, gross height in inches	10 11/16 (270)
External dimensions, gross width in inches	20 11/16 (525)
External dimensions, gross depth in inches	21 5/16 (540)
Net weight in lbs	10 (4.40)
Gross weight in lbs	15 (6.40)