


A 303  
Module  
for the optimum loading of up to 98 pipettes.



- Suitable for 98 measuring and full pipettes, for example
- Height of holding frame 150 mm
- 1 unit can be inserted into each of the baskets A 150 and APLW 106
- Vertical clearance w/out upper basket: 17 3/4" (450 mm) with A 150 lower basket

EAN: 4002515460825 / Material number: 09862560 / Old Material Number: 69630300D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PG 8583	•
PG 8583 CD	•
PG 8593	•
<b>Product category</b>	
Insert/module for pipettes	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black
<b>Application area</b>	
Reprocessing of pipettes/vials	•
<b>Capacity</b>	
Pipettes [number]	98
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in inches	7 5/16 (185)
External dimensions, net width in inches	8 3/4 (222)
External dimensions, net depth in inches	18 9/16 (471)
External dimensions, gross height in inches	10 11/16 (270)
External dimensions, gross width in inches	20 11/16 (525)
External dimensions, gross depth in inches	21 5/16 (540)
Net weight in lbs	6 (2.50)
Gross weight in lbs	8 (3.50)