

A 801  
 Adapter  
 for use of AUF 1 and AUF 2 in connection with A 105/1.



EAN: 4002515460931 / Material number: 09862940 /  
 Old Material Number: 69680100D

<b>Product category</b>	
Thermal disinfectors	•
<b>Machine affiliation</b>	
PG 8581	•
PG 8591	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
<b>Product category</b>	
Holder and/or adapter for transmission instruments	•
<b>Product characteristics</b>	
Material	Stainless steel
Color	Stainless Steel
<b>Application area</b>	
Reprocessing of dental instruments	•
<b>Dimensions and weight</b>	
External dimensions, net height in inches	2 (50)
External dimensions, net width in inches	2 (50)
External dimensions, net depth in inches	2 (50)
External dimensions, gross height in inches	0 1/2 (12)
External dimensions, gross width in inches	3 3/16 (80)
External dimensions, gross depth in inches	4 3/4 (120)
Net weight in lbs	1 (0.01)
Gross weight in lbs	1 (0.01)