

A 802

Jets

for the optimum cleaning of the dispenser compartment.



- For use in laboratory R+D machines with powder dispenser
- For use with injector modules
- Suitable for removing powder residues from dispenser compartment

EAN: 4002515460979 / Material number: 09863070 / Old Material Number: 69680200D

Product category	
Washer-disinfectors, laboratory	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8581	•
PG 8562	•
PG 8582	•
PG 8504	•
PG 8583	•
Product category	
Other accessories	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Dimensions and weight	
External dimensions, net height in inches	30
External dimensions, net width in inches	187
External dimensions, net depth in inches	15
External dimensions, gross height in inches	0 1/2 (12)
External dimensions, gross width in inches	3 3/16 (80)
External dimensions, gross depth in inches	4 3/4 (120)
Net weight in lbs	1 (0.02)
Gross weight in lbs	1 (0.02)