

A 802

Jets

for the optimum cleaning of the dispenser compartment.



- For use in laboratory R+D machines with powder dispenser
- For use with injector modules
- Suitable for removing powder residues from dispenser compartment

EAN: 4002515460979 / Material number: 09863070 / Old Material Number: 69680200D

Product category	
Washer-disinfectors, laboratory	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8581	•
PG 8562	•
PG 8582	•
PG 8504	•
PG 8583	•
Product category	
Other accessories	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Dimensions and weight	
External dimensions, net height in inches	1 3/16 (30)
External dimensions, net width in inches	7 3/8 (187)
External dimensions, net depth in inches	0 5/8 (15)
External dimensions, gross height in inches	0 1/2 (12)
External dimensions, gross width in inches	3 3/16 (80)
External dimensions, gross depth in inches	4 3/4 (120)
Net weight in lbs	1 (0.02)
Gross weight in lbs	1 (0.02)