

A 803
Holder
for the optimum loading of transmission instruments.



- For use in A 105/1 upper basket along with A 800 tubular filter
- Suitable for numerous contra-angles & hand pieces as well as turbines

EAN: 4002515460993 / Material number: 09863100 / Old Material Number: 69680300D

Product category	
Thermal disinfectors	•
Machine affiliation	
PG 8581	•
PG 8591	•
Product category	
Holder	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black
Application area	
Reprocessing of dental instruments	•
Dimensions and weight	
External dimensions, net height in inches	1 5/16 (32) in (mm)
External dimensions, net width in inches	1 5/16 (32) in (mm)
External dimensions, net depth in inches	1 5/8 (40) in (mm)
External dimensions, gross height in inches	3 3/16 (80)
External dimensions, gross width in inches	3 3/4 (95)
External dimensions, gross depth in inches	4 3/8 (110)
Net weight in lbs	1 (0.05)
Gross weight in lbs	1 (0.09)