

XKM 3000 L Med
Communication module
with Ethernet interface for integration in a network.



- Ring memory for max. 30 process reports and 2 graphic reports
- Additional components for direct network connection not required
- DHCP support for automatic setting of the network parameters
- Transfer speed 10/100 Mbit/s
- Length of power cord: 4.9 ft; cord extension possible

EAN: 4002515513392 / Material number: 09902230 / Old Material Number: 68822101D

Product category	
Washer-disinfectors, laboratory	•
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8583	•
PG 8583 CD	•
PG 8591	•
PG 8592	•
PG 8593	•
Product category	
Process documentation/operating documentation/communication	•
Product characteristics	
Material	Plastic
Color	Grey
Application area	
Reprocessing of pipettes/vials	•
Reprocessing of laboratory glassware	•
Reprocessing of slides	•
Reprocessing of butyrometers	•
Reprocessing of pharmacy/infusion bottles	•
Reprocessing of anesthetic instruments	•
Reprocessing of MIS instruments	•
Reprocessing of ophthalmic instruments	•
Reprocessing of robotic instruments	•
Reprocessing of GYN instruments	•
Reprocessing of ENT instruments	•
Reprocessing of dental instruments	•
Reprocessing of medical shoes	•
Reprocessing of basins	•
Reprocessing of kidney bowls	•
Reprocessing of baby bottles	•
Reprocessing of surgical instruments	•
Reprocessing of lumened instruments	•
Reprocessing of orthopedic-surgery instruments	•
Reprocessing of intramedullary nails/bone marrow drills	•
Reprocessing of trays	•
Reprocessing of mesh trays	•
Reprocessing of measuring cylinders	•
Reprocessing of containers, e.g., soda-lime containers	•
Reprocessing of different utensils	•
Dimensions and weight	
External dimensions, net height in inches	1 9/16 (39) in (mm)
External dimensions, net width in inches	3 1/2 (88) in (mm)
External dimensions, net depth in inches	4 3/8 (110) in (mm)

XKM 3000 L Med
Communication module
with Ethernet interface for integration in a network.



EAN: 4002515513392 / Material number: 09902230 / Old Material Number: 68822101D

External dimensions, gross height in inches	5 9/16 (140)
External dimensions, gross width in inches	7 1/2 (190)
External dimensions, gross depth in inches	11 5/8 (295)
Net weight in lbs	1 (0.44)
Gross weight in lbs	2 (0.61)
Standards, test certificates, and symbols	
CE	•