

A 800

Tubular filter

for the safe and economic reprocessing of hollow instruments.



- Filter in ply mesh stainless steel mesh for mobile units
- Mesh size 63 µm
- Reprocessable with cleaning brush A 804

EAN: 4002515538081 / Material number: 09960680 / Old Material Number: 68680000D

Product category	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
Product category	
Filter	•
Product characteristics	
Material	Stainless steel
Dimensions and weight	
External dimensions, net height in inches	0 15/16 (23) in (mm)
External dimensions, net width in inches	0 15/16 (23) in (mm)
External dimensions, net depth in inches	16 (405) in (mm)
External dimensions, gross height in inches	1 13/16 (45)
External dimensions, gross width in inches	8 11/16 (220)
External dimensions, gross depth in inches	17 9/16 (445)
Diameter in inches	0 15/16 (23)
Net weight in lbs	1 (0.36)
Gross weight in lbs	2 (0.88)