

DOS G 80/1

Hose dispenser pump f. 5 gal containers

With level indicator for the automatic dispensing of liquid cleaning agents.

EAN: 4002515529300 / Material number: 10001780 / Old Material Number: 67747025USA



- For freshwater dishwasher (PG8055-60) & tank dishwasher (G8066 & PG8067)
- 13 1/8 (333) in (mm) suction wand for 5 l media canisters
- Dispenses up to 55 Milliliter ml per minute
- Set via the washing machine controls

Product category	
Fresh-water dishwashers	•
Machine affiliation	
PG 8041	•
PG 8043	•
PG 8045	•
PG 8056	•
PG 8056 U	•
PG 8057 TD	•
PG 8057 TD U	•
PG 8058 U	•
PG 8059	•
PG 8059 U	•
PG 8061	•
PG 8061 U	•
PG 8096	•
PG 8096 U	•
PG 8099	•
PG 8099 U	•
Product category	
External dispenser system	•
Product characteristics	
Material	Plastic
Color	Grey
Production rate in ml/60 s	55
Length of suction wand in mm	13 1/8 (333)
Standard electrical connection	
Number of phases	1
Neutral conductor	N
AC/DC	AC
Voltage in V	200-240
Frequency in Hz	60
Total rated load in kW	0.01
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in inches	0
External dimensions, net width in inches	0
External dimensions, net depth in inches	0
External dimensions, gross height in inches	7 7/8 (200)
External dimensions, gross width in inches	18 15/16 (480)
External dimensions, gross depth in inches	5 15/16 (150)
Net weight in lbs	7 (2.80)
Gross weight in lbs	7 (3.00)