

## KWT 6312 UGS

Built-under wine storage unit

for perfect enjoyment and timeless design with its Push2open and SommelierSet.



- Perfect design for handleless kitchens – Push2open
- Odorless refrigeration section thanks to Active AirClean filters
- Present the storage of your wine in the perfect light – LED lighting
- Optimum wine protection through UV filters in tinted glass front
- Perfect enjoyment of your wine with the SommelierSet

EAN: 4002515539705 / Material number: 10028650 /  
Old Material Number: 36631200USA

General information	
Model number	KWT 6312 UGS
Manufacturer material number	10028650
EAN number	4002515539705
Class	Platinum
Appliance category	
Wine storage unit	•
Design	
Built-in	•
Built-in integrated	•
Built-in built-under	•
Door hinge side	right
Convertible hinging	•
Can be installed side-by-side	•
Design	
Front color	Glass door
Wine compartment lighting	LED illumination
User convenience	
DynaCool	•
SoftClose door closing damper	•
Active AirClean filter	•
Push2open	•
Total number of shelves	4
Number of shelves	4
Number of FlexiFrame Wire Racks	3
SommelierSet	•
NoteBoard	Yes
Low-vibration storage	•
Control	
Electronic temperature display and control	TouchControl
No. of temperature zones	1
Min. temperature setting range in °F (°C)	5
Safety	
Lock function	•
Acoustic door alarm	•
Acoustic temperature alarm	•
Optical door alarm	•
Technical data	
Niche width in inches (mm), min.	23 5/8 (600)
Niche width in inches (mm), max.	23 5/8 (600)
Niche height in inches (mm), min.	32 3/8 (820)
Niche height in inches (mm), max.	34 3/8 (870)
Niche depth in in.(mm)	22 7/8 (580)
Appliance width in in.(mm)	23 5/8 (597)
Appliance height in in.(mm)	32 3/8 (820)
Appliance depth in in.(mm)	22 3/4 (577)
Weight in lbs.(kg)	102 (46.0)
Wine compartment in L	3.88 (110)

## KWT 6312 UGS

Built-under wine storage unit

for perfect enjoyment and timeless design with its Push2open and SommelierSet.



EAN: 4002515539705 / Material number: 10028650 /


Old Material Number: 36631200USA

Total usable capacity in cu. ft. (l)	3.88 (110)
0.75 l Bordeaux bottle capacity (no. of bottles)	46 bottles
Storage time in case of fault in h	0
Current consumption in milliampere (mA)	1500
Voltage in V	115
Fuse rating in A	10
Number of phases	1
Frequency in Hz	60
Wire length in ft (m)	7 (2.4)

KWT 6312 UGS

Built-under wine storage unit

for perfect enjoyment and timeless design with its Push2open and SommelierSet.



KWT6312UGS, Interior dimensions (Installation drawings)