


UG 30-90/60-85

Box toe-kick

to raise the machine incl. side unit, 23.6 inches deep, 35.4 inches wide.



- Usable for Aqua Purificator PG 8595 and DOS cabinet PG 8596

EAN: 4002515558348 / Material number: 10087360 / Old Material Number: 69309001D

Product category	
Washer-disinfectors, laboratory	•
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8581	•
PG 8562	•
PG 8582	•
PG 8504	•
PG 8583	•
PG 8591	•
PG 8592	•
PG 8593	•
Product category	
Plinth/stand	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Dimensions and weight	
External dimensions, net height in inches	11 13/16 (300)
External dimensions, net width in inches	35 3/8 (897)
External dimensions, net depth in inches	21 3/8 (542)
External dimensions, gross height in inches	13 (330)
External dimensions, gross width in inches	36 1/16 (915)
External dimensions, gross depth in inches	22 3/16 (563)
Net weight in lbs	30 (13.40)
Gross weight in lbs	34 (15.00)