

UBS 1

Conversion kit

For use with selected series G 78 load carriers in PG 85 machines.



- For older loading trolleys with dryer connection

EAN: 4002515588970 / Material number: 10131090 / Old Material Number: 69215014D

Product category	
Washer-disinfectors, laboratory	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8583	•
PG 8583 CD	•
PG 8592	•
PG 8593	•
Product category	
Conversion kit	•
Product characteristics	
Material	Plastic
Color	Black
Dimensions and weight	
External dimensions, net height in inches	3 (75)
External dimensions, net width in inches	9 5/16 (235)
External dimensions, net depth in inches	6 1/2 (165)
External dimensions, gross height in inches	3 3/16 (80)
External dimensions, gross width in inches	9 1/2 (240)
External dimensions, gross depth in inches	6 3/4 (170)
Net weight in lbs	1 (0.10)
Gross weight in lbs	1 (0.17)