

MC/2

Transfer trolley

height adjustable, for loading & unloading washer-&disinfectors.



- Equipped with 4 castors, 2 of which can be locked
- 2 loading racks (trough shaped) as well as handle and cover plate

EAN: 4002515588628 / Material number: 10138690 / Old Material Number: 69200109D

<b>Product category</b>	
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8535	•
PG 8536	•
PG 8581	•
PG 8582	•
PG 8582 CD	•
PG 8583	•
PG 8583 CD	•
PG 8591	•
PG 8592	•
PG 8593	•
<b>Product category</b>	
Storage/transport system	•
<b>Product characteristics</b>	
Material	Stainless steel
Color	Stainless Steel
<b>Dimensions and weight</b>	
External dimensions, net height in inches	39 1/2 (1003)
External dimensions, net width in inches	23 5/8 (600)
External dimensions, net depth in inches	31 3/4 (805)
External dimensions, gross height in inches	44 1/16 (1119)
External dimensions, gross width in inches	32 7/8 (835)
External dimensions, gross depth in inches	50 7/16 (1280)
Net weight in lbs	76 (34.20)
Gross weight in lbs	116 (52.40)