

GP CL KM 0252 L

Cer. and stainless cleaner 8.5 fl oz.

For best cleaning results and safe use.

EAN: 4002515606056 / Material number: 10173130 / Old Material Number: 22996291EU2



- Best product on the market –Test winner in Switzerland
- Removes grease, baked-on residues and water marks powerfully
- Leaves a protective coating, making cleaning easier the next time
- Gentle cleaning thanks to fine-crystal formula
- Do not use with Clean Touch Steel™ or PerfectClean
- Optimum care for many years of reliability