


UG 30-90/70-85

Box toe-kick

to raise the machine, 27.6 inches deep, 35.4 inches wide.



- Suitable for all machines, width 35.4, depth 27.6 inches

EAN: 4002515622506 / Material number: 10210600 / Old Material Number: 69309002D

Product category	
Washer-disinfectors, laboratory	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8582 CD	•
PG 8583 CD	•
Product category	
Plinth/stand	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Dimensions and weight	
External dimensions, net height in inches	12 1/2 (317)
External dimensions, net width in inches	35 7/16 (900)
External dimensions, net depth in inches	27 9/16 (700)
External dimensions, gross height in inches	12 1/4 (310)
External dimensions, gross width in inches	24 1/16 (610)
External dimensions, gross depth in inches	24 1/16 (610)
Net weight in lbs	35 (15.60)
Gross weight in lbs	38 (17.00)