

UGT 44-60/85  
 Box toe-kick, with door  
 for ergonomic loading and unloading of the machine.



- Offers the option of storing e.g., process chemicals

EAN: 4002515647295 / Material number: 10259260 / Old Material Number: 69506001D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Thermal disinfectors	•
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8581	•
PG 8582	•
PG 8583	•
PG 8591	•
PG 8592	•
PG 8593	•
<b>Product category</b>	
Plinth/stand	•
<b>Product characteristics</b>	
Material	Stainless steel
Color	Stainless Steel
<b>Dimensions and weight</b>	
External dimensions, net height in inches	17 3/8 (440)
External dimensions, net width in inches	23 5/8 (600)
External dimensions, net depth in inches	23 5/8 (600)
External dimensions, gross height in inches	18 1/8 (460)
External dimensions, gross width in inches	24 1/16 (610)
External dimensions, gross depth in inches	24 1/16 (610)
Net weight in lbs	12 (5.00)
Gross weight in lbs	15 (6.40)