

ProCare Lab 10 AP - 10 l

Liquid detergent, alkaline, 10 l

for optimum reprocessing of laboratory utensils.



- Free from surfactants, phosphorus compounds, and oxidation agents
- Specially developed for Miele Professional
- Suitable for all types of materials
- pH value 12.3-12.5

EAN: 4002515660713 / Material number: 10265640 / Old Material Number: 62924201EU1

Product category	
Washer-disinfectors, laboratory	•
Applications	
Removal of low-viscosity oils	Good
Removal of inorganic and organic soiling	Good
Removal of stubborn or dried-on residues	Good
Dispensing recommendation	
Dispensing recommendation per rinse	2-4 ml/l (0.2-0.4%)
pH value of solution	12.3-12.5
Main applications	
Medical laboratories	•
Biological laboratories	•
Microbiological laboratories	•
Chemical laboratories	•
Water laboratories	•
Laboratories for the food industry	•
Laboratories for the cosmetics industry	•
Laboratories for the mineral-oil industry	•
Phosphate industry	•
Laboratories for the pharmaceutical industry	•
Product characteristics	
Medium type	Alkaline cleaning agent
Aggregate form	Liquid
Ingredients	Alkali hydroxide, Polycarboxylate
Suitable for plastic	•
Suitable for stainless steel	•
Suitable for laboratory glassware	•
Suitable for ceramics	•
Not suitable for anodized aluminum	•
Does not contain oxidants	•
Does not contain phosphates	•
Does not contain surfactants	•
Oxidizing effect	•
Not suitable for light alloys	•
Not suitable for non-anodized aluminum	•
Not suitable for light and non-ferrous metals	•
Dimensions and weight	
External dimensions, net height in inches	12 3/16 (308)
External dimensions, net width in inches	7 9/16 (191)
External dimensions, net depth in inches	9 1/4 (234)
External dimensions, gross height in inches	12 5/8 (320)
External dimensions, gross width in inches	9 1/16 (230)
External dimensions, gross depth in inches	7 1/2 (190)
Gross weight in lbs	0.49 (14)
Fill volume in ml	338.15 (10000) fl oz (US) (Milliliter)
Compliance, certificates/symb.	
GHS/CLP compliant	•