

ProCare Lab 10 MA - 5 l

Liquid detergent, mildly alkaline, 5 l

for optimum reprocessing of laboratory glassware and utensils.

EAN: 4002515660812 / Material number: 10265670 / Old Material Number: 62923101EU1



- Free from surfactants and oxidants
- Specially developed for Miele Professional
- Suitable for all types of materials
- PH value 11.6-11.8

Product category	
Washer-disinfectors, laboratory	•
Applications	
Removal of blood	Good
Removal of protein	Good
Removal of residues of nutrient mediums	Good
Removal of organic residues	Good
Dispensing recommendation	
Dispensing recommendation per rinse	2-4 ml/l (0.2-0.4%)
pH value of solution	11.6-11.8
Main applications	
Medical laboratories	•
Biological laboratories	•
Chemical laboratories	•
Product characteristics	
Medium type	Alkaline cleaning agent
Aggregate form	Liquid concentrate
Ingredients	Alkali
Suitable for plastic	•
Suitable for stainless steel	•
Suitable for laboratory glassware	•
Suitable for ceramics	•
Limited suitability for non-ferrous metals	•
Not suitable for anodized aluminum	•
Does not contain oxidants	•
Does not contain surfactants	•
Limited suitability for non-anodized aluminum	•
Dimensions and weight	
External dimensions, net height in inches	8 9/16 (217)
External dimensions, net width in inches	5 5/8 (142)
External dimensions, net depth in inches	9 13/16 (248)
External dimensions, gross height in inches	9 7/8 (250)
External dimensions, gross width in inches	7 1/2 (190)
External dimensions, gross depth in inches	5 15/16 (150)
Gross weight in lbs	0.18 (5)
Fill volume in ml	169.07 (5000) fl oz (US) (Milliliter)
Compliance, certificates/symb.	
GHS/CLP compliant	•