

ProCare Lab 11 AP - 10 kg  
 Powder cleaning agent, alkaline, 10 kg  
 for optimum reprocessing of laboratory utensils.

EAN: 4002515650868 / Material number: 10266540 / Old Material Number: 62922101EU1



- Free from phosphorous compounds, oxidising effect
- Specially developed for Miele Professional
- Suitable for all types of materials
- PH value 11.5-11.7

<b>Product category</b>	
Washer-disinfectors, laboratory	•
<b>Applications</b>	
Removal of inorganic and organic soiling	Good
<b>Dispensing recommendation</b>	
pH value of solution	11.5-11.7
Recommended water range in mmol/l	0-1.1
<b>Main applications</b>	
Water laboratories	•
Serological institutes	•
Phosphate industry	•
<b>Product characteristics</b>	
Medium type	Alkaline cleaning agent
Aggregate form	Powder
Ingredients	Oxygen bleach, Tensides
Suitable for plastic	•
Suitable for stainless steel	•
Suitable for porcelain/ceramics	•
Suitable for laboratory glassware	•
Not suitable for anodized aluminum	•
Does not contain phosphates	•
Not suitable for non-anodized aluminum	•
Not suitable for light and non-ferrous metals	•
<b>Dimensions and weight</b>	
External dimensions, net height in inches	10 9/16 (268)
External dimensions, net width in inches	10 7/16 (265)
External dimensions, net depth in inches	10 7/16 (265)
External dimensions, gross height in inches	11 1/16 (280)
External dimensions, gross width in inches	10 11/16 (270)
External dimensions, gross depth in inches	10 11/16 (270)
Load size in g	352.73 (10000)
Gross weight in lbs	0.35 (10)
<b>Compliance, certificates/symb.</b>	
GHS/CLP compliant	•