

PPSSW0010

Thermoelectric test

for washing machines and dishwashers.

- Record temperatures and holding times
- Use of calibrated data loggers
- Test measurement protocols

EAN: 4002515756737 / Material number: 10486590 / Old Material Number:

<b>Product category</b>	
Washing machines	•
Fresh-water dishwashers	•
<b>Contract contents</b>	
Obtaining the necessary documentation	•
Traveling for selected service	•
Working time for selected service	•
Detailed documentation	•
Thermoelectric test using calibrated measuring equipment	•
Interval	1 x jährlich