

## A 815 Nozzle

For perfect reprocessing of transmission instruments.

EAN: 4002515830413 / Material number: 10619690 / Old Material Number: 68681500D



- For use in A 812, A 813, or in A 105/1 upper basket
- Consists of 4 connectors with silicone hoses, L approx. 4,5"

<b>Product category</b>	
Thermal disinfectors	•
<b>Machine affiliation</b>	
PG 8581	•
PG 8591	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
<b>Product category</b>	
Injector nozzle	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Color	Stainless Steel, Transparent
<b>Application area</b>	
Reprocessing of dental instruments	•
<b>Capacity</b>	
Transmission instruments [number]	4
<b>Dimensions and weight</b>	
External dimensions, net height in inches	0 7/16 (10)
External dimensions, net width in inches	0 7/16 (10)
External dimensions, net depth in inches	4 1/2 (114)
External dimensions, gross height in inches	3 3/16 (80)
External dimensions, gross width in inches	3 3/4 (95)
External dimensions, gross depth in inches	4 3/8 (110)
Diameter in inches	0 7/16 (10)
Net weight in lbs	1 (0.05)
Gross weight in lbs	1 (0.12)