

A 833

Injector nozzle

for optimum cleaning of a variety of wash items.



- Suitable for mobile injector units, with screw connection
- For optimum cleaning results thanks to lateral drill holes

EAN: 4002515941225 / Material number: 10635480 / Old Material Number: 69745302D

<b>Product category</b>	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
G 7831	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
PWA 8672	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
<b>Product category</b>	
Injector nozzle	•
<b>Product characteristics</b>	
Material	Stainless steel
Color	Stainless Steel
<b>Application area</b>	
Reprocessing of MIS instruments	•
Reprocessing of dental instruments	•
Reprocessing of lumened instruments	•
<b>Dimensions and weight</b>	
External dimensions, net width in inches	0 3/16 (4)
External dimensions, net depth in inches	4 3/8 (110)
External dimensions, gross height in inches	0 1/2 (12)
External dimensions, gross width in inches	3 3/16 (80)
External dimensions, gross depth in inches	4 3/4 (120)
Net weight in lbs	1 (0.00)
Gross weight in lbs	1 (0.01)