

A 816

Injector nozzle

for optimum cleaning of a variety of wash items.



- Suitable for mobile injector units and baskets, screw connection
- For optimum cleaning results thanks to lateral drill holes

EAN: 4002515978849 / Material number: 10635660 /
Old Material Number: 69681600D

Product category	
Thermal disinfectors	•
Cleaning machine, small, dental	•
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
G 7831	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
PWA 8672	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Injector nozzle	•
Application area	
Reprocessing of MIS instruments	•
Reprocessing of dental instruments	•
Reprocessing of lumened instruments	•
Product characteristics	
Material	Stainless steel
Color	Stainless steel
Dimensions and weight	
External dimensions, net width in inches	0 1/8 (3)
External dimensions, net depth in inches	3 9/16 (90)
External dimensions, gross height in inches	0 1/2 (12)
External dimensions, gross width in inches	3 3/16 (80)
External dimensions, gross depth in inches	4 3/4 (120)
Net weight in lbs	0 (0.01)
Gross weight in lbs	0 (0.01)