

Classic C1 Cat & Dog PowerLine - SBBN0
 canister vacuum cleaners
 with electrobrush for thorough cleaning of heavy-duty carpeting.



- High suction power – 1200 W
- Optimum carpet care – electrobrush for intensive deep cleaning
- Adjustable working height – Stainless steel telescopic tube
- Effortless vacuuming of large areas w/ 29.5 ft operating radius
- Odors are neutralized more effectively – Active AirClean filter

EAN: 4002515818091 / Material number: 10639470 / Old Material Number: 41BBN031USA

Design	
Canister vacuum cleaner with dustbag	•
Product line	Classic C1
Motor type	PowerLine
Product type	SBBN0
Construction type and design	
Color	Lotus white
User convenience	
Electronic suction power setting	6-stage rotary selector
Electrobrush socket	•
Ergonomic carrying handle	•
Ergonomic handle	•
Comfort handle	•
Silence system	•
Telescopic suction wand	Electro telescopic wand
Park system for vacuuming breaks	•
Park system for storage	on one side
Mobility	
Caster wheels gentle on floors	•
Three-point maneuverability	•
Operating radius in ft (m)	29 (9)
Cord length in ft (m)	18 (5.5)
Filter system	
FilterBag type	HyClean GN
FilterBag™ volume in qt (US) (l)	4.8 (4.5) qt (US) (l)
FilterBag™ change indicator	•
Exhaust filter	AirClean filter
Motor protection filter	•
Safety	
Locking system for floorhead and handle	•
Automatic FilterBag™ positioning	•
Technical data	
Power	1200
Width in inches (mm)	18 1/4 (463)
Height in inch (mm)	279
Weight in lbs (kg)	21 (9.1)
Power rating in W	1000
Standard accessories	
Three-piece accessory set	On VarioClip
Upholstery nozzle	•
Crevice nozzle	•
Dusting brush with synthetic bristles	•
Floorbrush	SBB Parquet-3 Small Parquet Brush
Small Parquet Floorhead SBB Parquet-3	•
Hand turbobrush	STB 101
Handheld Mini Turbo Brush STB 101	•
Electrobrush	SEB 228 Electro Plus
Electro Plus Floorhead SEB 228	•
Filter	AirClean filter
Cleaning brush	No