

DTD 2
Drip tray
for the secure placement of ProCare media canisters.

EAN: 4002515933695 / Material number: 10695460 / Old Material Number: 61910001D

Dimensions and weight	
External dimensions, gross height in inches	2 15/16 (74)
External dimensions, gross width in inches	12 7/8 (327)
External dimensions, gross depth in inches	18 1/4 (463)
Gross weight in lbs	0.04 (1)



- Space for 2 ProCare media canisters
- Secure place for 5 and 10 l canisters in the ProCare series
- Optimum residual emptying of the canister due to angle
- Suitable for use in an adjacent cabinet