

A 819

Irrigation sleeves

for optimum cleaning of lumened instruments.

EAN: 4002515990278 / Material number: 10721750 / Old Material Number: 68681900D



- E.g., for magnetostrictive ultrasound scaler tips, e.g., Cavitron
- Int. dimension 4 1/4" (108 mm)
- Including cap with an opening of Ø 3/8" (10 mm)
- With mounting clip

| | |
|---|--------------------------|
| Product category | |
| Thermal disinfectors | • |
| Machine affiliation | |
| PG 8528 | • |
| G 7831 | • |
| PG 8536 | • |
| PG 8581 | • |
| PG 8591 | • |
| Product category | |
| Holder | • |
| Product characteristics | |
| Material | Stainless steel, Plastic |
| Color | Stainless Steel, Black |
| Application area | |
| Reprocessing of dental instruments | • |
| Capacity | |
| Ultrasonic tips [number] | 1 |
| Dimensions and weight | |
| External dimensions, net height in inches | 0 7/16 (10) in (mm) |
| External dimensions, net width in inches | 3 15/16 (100) in (mm) |
| External dimensions, net depth in inches | 4 1/2 (114) in (mm) |
| External dimensions, gross height in inches | 0 1/2 (12) |
| External dimensions, gross width in inches | 4 3/4 (120) |
| External dimensions, gross depth in inches | 6 3/4 (170) |
| Net weight in lbs | 1 (0.02) |
| Gross weight in lbs | 1 (0.02) |