

DOS G 80 Flex
Hose metering pump
with level indicator for the automatic dispensing of liquid cleaning agents.



- For canisters up to 13 3/8" (340 mm) tall (to canister opening)
- Dispenses up to 55 Milliliter ml per minute
- Set via the dishwasher controls
- Flexible vent valve cap
- 13 7/8 (352) in (mm) inch suction wand for 5 and 10 l canisters

EAN: 4002516280064 / Material number: 10730850 / Old Material Number: 67747040USA

Product category	
Fresh-water dishwashers	•
Tank dishwashers	•
Washer-disinfectors, laboratory	•
Machine affiliation	
PG 8056	•
PG 8056 U	•
PG 8059	•
PG 8059 U	•
PG 8060	•
PG 8061	•
PG 8061 U	•
PG 8583	•
PG 8593	•
Product category	
External dispenser system	•
Product characteristics	
Material	Plastic
Color	Grey
Production rate in ml/60 s	55
Length of suction wand in mm	13 7/8 (352)
Length of dispensing hose to suction wand in inches (mm)	1800
Length of dispensing hose from machine to dispenser module in inches (mm)	2850
Length of power cord to machine in inches (mm)	2800
Minimum width of canister opening in inches (mm)	31
Maximum width of canister opening in inches (mm)	60
Standard electrical connection	
Number of phases	1
Neutral conductor	N
AC/DC	AC
Voltage in V	200-240
Frequency in Hz	60
Total rated load in kW	0.01
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in inches	4 1/2 (114) in (mm)
External dimensions, net width in inches	4 9/16 (115) in (mm)
External dimensions, net depth in inches	5 1/2 (139) in (mm)
External dimensions, gross height in inches	6 5/16 (160)
External dimensions, gross width in inches	7 13/16 (198)
External dimensions, gross depth in inches	18 7/8 (478)
Net weight in lbs	6 (2.27)
Gross weight in lbs	6 (2.51)