


A 105/1

Upper basket

Height adjustable, without injector nozzles, for various instruments.



- Right-hand side: space for 22 injectors
- Left-hand side free for inserts
- Active inner drying of hollow instruments
- Active inner drying of dental and surgical handpieces
- A 800 tubular filter can be used in both injector manifolds

EAN: 4002515941157 / Material number: 10777260 / Old Material Number: 69610501D

Product category	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
Product category	
Injector basket	•
Upper basket	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black
Option to use a central filter	•
Application area	
Reprocessing of ENT instruments	•
Reprocessing of dental instruments	•
Capacity	
Transmission instruments [number]	22
Lumened instruments [number]	22
Dimensions and weight	
External dimensions, net height in inches	8 (202)
External dimensions, net width in inches	20 13/16 (528)
External dimensions, net depth in inches	20 13/16 (528)
External dimensions, gross height in inches	10 11/16 (270)
External dimensions, gross width in inches	20 11/16 (525)
External dimensions, gross depth in inches	21 5/16 (540)
Net weight in lbs	11 (4.86)
Gross weight in lbs	13 (5.56)
Accessory pack	
A 838 key set [number]	1