

TC 900

Dresser panelling

For machines 36 in. (900 mm) wide without a steam condenser or drying unit.



- Includes service door with lock
- Diaphragm wall paneling up to ceiling must be provided on site

EAN: 4002515991756 / Material number: 10880700 / Old Material Number: 69252003D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8615	•
PLW 8616	•
Product category	
Cover	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Dimensions and weight	
External dimensions, net height in inches	33 1/8 (840)
External dimensions, net width in inches	33 5/16 (845)
External dimensions, net depth in inches	41 3/16 (1045)
External dimensions, gross height in inches	39 (990)
External dimensions, gross width in inches	33 5/16 (845)
External dimensions, gross depth in inches	47 1/4 (1200)
Net weight in lbs	108 (48.90)
Gross weight in lbs	141 (63.90)