

TC 1150

Dresser panelling

For machines 45 in. (1150 mm) wide without a steam condenser or drying unit.

EAN: 4002515992128 / Material number: 10880750 / Old Material Number: 69252004D



- Includes service door with lock
- Diaphragm wall paneling up to ceiling must be provided on site

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8617	•
Product category	
Cover	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Dimensions and weight	
External dimensions, net height in inches	39 1/2 (1002)
External dimensions, net width in inches	32 7/8 (835)
External dimensions, net depth in inches	47 1/4 (1200)
External dimensions, gross height in inches	39 1/2 (1002)
External dimensions, gross width in inches	32 7/8 (835)
External dimensions, gross depth in inches	47 1/4 (1200)
Net weight in lbs	133 (60.13)
Gross weight in lbs	166 (75.13)