

DS 10

Diaphragm dispensing pump

For 10-liter media canisters, for dispensing liquid media.



- For retrofitting additional dispenser pumps
- For dispensing detergents, rinse aid, and neutralizing agents
- Includes a suction wand for 10-liter media canister
- Scope of delivery: Dispenser pump, lance, volume control, & hoses

EAN: 4002515992241 / Material number: 10880870 / Old Material Number: 69249101D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Dispenser system	•
Standard electrical connection	
Number of phases	0
AC/DC	DC
Voltage in V	24
Frequency in Hz	0
Total rated load in kW	0.01
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in inches	0
External dimensions, net width in inches	0
External dimensions, net depth in inches	0
External dimensions, gross height in inches	7 7/8 (200)
External dimensions, gross width in inches	19 (480)
External dimensions, gross depth in inches	6 (150)
Net weight in lbs	3 (1.05)
Gross weight in lbs	3 (1.25)