

A 501

Mobile unit

With 4 levels for cleaning laboratory glassware and utensils.



- For holding inserts and mesh trays
- Extendable to 5 or 6 levels with A 601 racks

EAN: 4002516050438 / Material number: 10973890 / Old Material Number: 69650101D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 8615	•
PLW 8616	•
PLW 8617	•
<b>Product category</b>	
Mobile unit for inserts	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
DIN mesh trays [number]	18
SPRI mesh trays [number]	12
<b>Dimensions and weight</b>	
External dimensions, net height in inches	26 (659)
External dimensions, net width in inches	25 5/8 (650)
External dimensions, net depth in inches	28 (710)
External dimensions, gross height in inches	31 5/16 (795)
External dimensions, gross width in inches	33 7/8 (860)
External dimensions, gross depth in inches	28 3/16 (715)
Net weight in lbs	22 (9.80)
Gross weight in lbs	69 (31.00)