

A 503

Mobile unit

No levels, for cleaning laboratory glassware/utensils.



- For loading modules and inserts
- Depending on the height, up to 6 modules can be used
- Upper rail and upper strut can be removed

EAN: 4002516050445 / Material number: 10973900 /
Old Material Number: 69650301D

Product category

Large-capacity washer-disinfectors, laboratory •

Machine affiliation

PLW 8615 •

PLW 8616 •

PLW 8617 •

Product category

Mobile unit with module connection •

Product characteristics

Material Stainless steel, Plastic

Color Stainless Steel, Black

Application area

Reprocessing of laboratory glassware •

Dimensions and weight

External dimensions, net height in inches 17 1/16 (433)

External dimensions, net width in inches 31 1/8 (790)

External dimensions, net depth in inches 25 5/8 (650)

External dimensions, gross height in inches 31 1/8 (790)

External dimensions, gross width in inches 33 7/8 (860)

External dimensions, gross depth in inches 28 3/16 (715)

Net weight in lbs 38 (17.00)

Gross weight in lbs 45 (20.00)