

A 503

Mobile unit

No levels, for cleaning laboratory glassware/utensils.



- For holding inserts and mesh trays
- Depending on the height, up to 6 modules can be used
- Upper rail and upper strut can be removed

EAN: 4002516050445 / Material number: 10973900 / Old Material Number: 69650301D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Mobile unit with module connection	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in inches	17 1/16 (433)
External dimensions, net width in inches	31 1/8 (790)
External dimensions, net depth in inches	25 5/8 (650)
External dimensions, gross height in inches	31 1/8 (790)
External dimensions, gross width in inches	33 7/8 (860)
External dimensions, gross depth in inches	28 3/16 (715)
Net weight in lbs	38 (17.00)
Gross weight in lbs	45 (20.00)