

A 606  
Module

For use of e.g., lab glassware inserts on the top level of A 503 mobile units



- With spray arm

EAN: 4002516050537 / Material number: 10973990 / Old Material Number: 69660601D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 8615	•
PLW 8616	•
PLW 8617	•
<b>Product category</b>	
Other inserts/other modules	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Dimensions and weight</b>	
External dimensions, net height in inches	3 7/8 (97)
External dimensions, net width in inches	19 15/16 (505)
External dimensions, net depth in inches	23 7/8 (606)
External dimensions, gross height in inches	4 3/4 (120)
External dimensions, gross width in inches	20 11/16 (525)
External dimensions, gross depth in inches	24 7/16 (620)
Net weight in lbs	6 (2.40)
Gross weight in lbs	7 (3.00)