


A 613  
Module

For the optimum reprocessing of volumetric and measuring pipettes.



- Equipped with 98 injector nozzles
- Suitable for measuring pipettes with a volume of 1–50 ml
- Suitable for volumetric pipettes with a volume of 1–100 ml
- For use in A 503 mobile unit

EAN: 4002516063360 / Material number: 10999070 / Old Material Number: 69661300D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 8615	•
PLW 8616	•
PLW 8617	•
<b>Product category</b>	
Insert/module for pipettes	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black
<b>Application area</b>	
Reprocessing of pipettes/vials	•
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Pipettes [number]	98
<b>Dimensions and weight</b>	
External dimensions, net height in inches	7 11/16 (194)
External dimensions, net width in inches	10 (254)
External dimensions, net depth in inches	23 13/16 (604)
External dimensions, gross height in inches	7 7/8 (200)
External dimensions, gross width in inches	10 11/16 (270)
External dimensions, gross depth in inches	24 1/16 (610)
Net weight in lbs	10 (4.40)
Gross weight in lbs	15 (6.40)