

A 620  
Module

For the optimum reprocessing of laboratory glassware, 10 injector positions.



- Can be equipped with up to 10 nozzles based on individual needs
- For use in A 503 mobile unit
- Part of the SimpleLoad system

EAN: 4002516063391 / Material number: 10999100 / Old Material Number: 69662000D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 8615	•
PLW 8616	•
PLW 8617	•
<b>Product category</b>	
Insert/module for narrow-neck glassware	•
Insert for various bottles	•
<b>Product characteristics</b>	
Material	Stainless steel
Color	Stainless Steel
Part of the SimpleLoad system	•
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Dimensions and weight</b>	
External dimensions, net height in inches	4 1/2 (113)
External dimensions, net width in inches	5 1/4 (133)
External dimensions, net depth in inches	23 15/16 (608)
External dimensions, gross height in inches	5 3/4 (145)
External dimensions, gross width in inches	5 15/16 (150)
External dimensions, gross depth in inches	24 7/16 (620)
Net weight in lbs	2 (0.77)
Gross weight in lbs	5 (1.82)