

A 622
Module

For the optimum reprocessing of laboratory glassware, 36 injector positions.



- Can be equipped with up to 36 nozzles based on individual needs
- For use in A 503 mobile unit
- Part of the SimpleLoad system

EAN: 4002516063414 / Material number: 10999120 / Old Material Number: 69662200D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Insert/module for narrow-neck glassware	•
Insert for various bottles	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Part of the SimpleLoad system	•
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Narrow-neck glassware [number]	36
Flasks [number]	36
Round flasks [number]	36
Dimensions and weight	
External dimensions, net height in inches	4 1/2 (113)
External dimensions, net width in inches	8 1/16 (204)
External dimensions, net depth in inches	23 15/16 (608)
External dimensions, gross height in inches	5 3/4 (145)
External dimensions, gross width in inches	8 11/16 (220)
External dimensions, gross depth in inches	24 5/8 (625)
Net weight in lbs	4 (1.42)
Gross weight in lbs	6 (2.53)