

A 846

Injector nozzle

With support, for cleaning laboratory flasks and round flasks.



- For laboratory bottles up to 10 L and round flasks up to 6 L
- Nozzle length: 5" (127 mm), nozzle length with holder: 6 1/4" (158 mm)
- Nozzle Ø: 10 mm
- Suitable for A 313 and A 612 modules

EAN: 4002516071464 / Material number: 11025590 / Old Material Number: 69684600D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
<b>Product category</b>	
Injector nozzle	•
<b>Product characteristics</b>	
Material	Stainless steel, null, Plastic
Color	Stainless Steel, Black, Grey
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Dimensions and weight</b>	
External dimensions, net height in inches	6 1/4 (158) in (mm)
External dimensions, net width in inches	4 3/8 (110) in (mm)
External dimensions, net depth in inches	4 3/8 (110) in (mm)
External dimensions, gross height in inches	6 5/16 (160)
External dimensions, gross width in inches	5 3/4 (145)
External dimensions, gross depth in inches	10 7/16 (265)
Net weight in lbs	1 (0.23)
Gross weight in lbs	1 (0.37)