

A 848

Injector nozzle

With support, for cleaning measuring flasks.



- For measuring flasks, up to 5 L
- Nozzle length: 1' 3/4" (350 mm)
- Nozzle Ø: 3/8" (10 mm), Ø of support: 5 9/16" (142 mm)
- Suitable for A 313 and A 612 modules

EAN: 4002516071488 / Material number: 11025610 / Old Material Number: 69684800D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Injector nozzle	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Stainless Steel, Grey
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in inches	13 13/16 (350) in (mm)
External dimensions, net width in inches	0 13/16 (20) in (mm)
External dimensions, net depth in inches	0 13/16 (20) in (mm)
External dimensions, gross height in inches	13 13/16 (350)
External dimensions, gross width in inches	7 1/8 (180)
External dimensions, gross depth in inches	2 (50)
Net weight in lbs	1 (0.30)
Gross weight in lbs	2 (0.53)