

A 306/1 Module

For the optimum loading of large-volume laboratory glassware.



- Suitable for 4 1–2 liter measuring cylinders, high shape
- Suitable for 2 1–2 liter measuring cyl., low shape + 2 high shape
- Plastic-coated contact surfaces
- Measuring cylinders even removable when module is in machine
- Maximum length of the measuring cylinders in the A 150 lower basket: 500 mm
- Up to 2 modules can be inserted in baskets A 150 and APL 106

EAN: 4002516076254 / Material number: 11035210 / Old Material Number: 69630601D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 7111	•
Product category	
Insert/module for measuring cylinder	•
Module for large-volume laboratory glassware	•
Application area	
Reprocessing of laboratory glassware	•
Reprocessing of measuring cylinders	•
Product characteristics	
Material	Stainless steel, Plastic, Wire, rilsan-coated
Color	Stainless Steel, Black, White
Capacity	
Measuring cylinders [number]	4
Dimensions and weight	
External dimensions, net height in inches	9 5/16 (235)
External dimensions, net width in inches	16 1/2 (418)
External dimensions, net depth in inches	18 9/16 (471)
External dimensions, gross height in inches	10 11/16 (270)
External dimensions, gross width in inches	20 11/16 (525)
External dimensions, gross depth in inches	21 5/16 (540)
Net weight in lbs	5 (2.10)
Gross weight in lbs	7 (3.10)