

A 306/1

Module

For the optimum loading of large-volume laboratory glassware.



- Suitable for 4 1–2 liter measuring cylinders, high shape
- Suitable for 2 1–2 liter measuring cyl., low shape + 2 high shape
- Plastic-coated contact surfaces
- Measuring cylinders even removable when module is in machine
- Maximum length of measuring cylinders: 19 11/16" (500 mm)
- Up to 2 modules can be inserted in baskets A 150 and APL 106

EAN: 4002516076254 / Material number: 11035210 / Old Material Number: 69630601D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
Product category	
Insert/module for measuring cylinder	•
Module for large-volume laboratory glassware	•
Product characteristics	
Material	Stainless steel, Plastic, Wire, rilsan-coated
Color	Stainless Steel, Black, White
Application area	
Reprocessing of laboratory glassware	•
Reprocessing of measuring cylinders	•
Capacity	
Measuring cylinders [number]	4
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in inches	235
External dimensions, net width in inches	418
External dimensions, net depth in inches	471
External dimensions, gross height in inches	10 11/16 (270)
External dimensions, gross width in inches	20 11/16 (525)
External dimensions, gross depth in inches	21 5/16 (540)
Net weight in lbs	5 (2.10)
Gross weight in lbs	7 (3.10)