

APH 592
 Colour transfer belt
 for protocol printer APH 510 PRT110.



EAN: 8809166675907 / Material number: 11053500 / Old Material Number: 68817005D

Product category	
Washer-disinfectors, laboratory	•
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
G 7831	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8583	•
PG 8583 CD	•
PG 8591	•
PG 8592	•
PG 8593	•
PWD 8541	•
PWD 8545	•
PWD 8546	•
PWD 8549	•
Product category	
Process documentation/operating documentation/ communication	•
Product characteristics	
Color	Red, Black
Application area	
Reprocessing of pipettes/vials	•
Reprocessing of laboratory glassware	•
Reprocessing of slides	•
Reprocessing of butyrometers	•
Reprocessing of pharmacy/infusion bottles	•
Reprocessing of anesthetic instruments	•
Reprocessing of MIS instruments	•
Reprocessing of ophthalmic instruments	•
Reprocessing of robotic instruments	•
Reprocessing of GYN instruments	•
Reprocessing of ENT instruments	•
Reprocessing of dental instruments	•
Reprocessing of medical shoes	•
Reprocessing of basins	•
Reprocessing of kidney bowls	•
Reprocessing of baby bottles	•
Reprocessing of surgical instruments	•
Reprocessing of lumened instruments	•
Reprocessing of orthopedic-surgery instruments	•
Reprocessing of intramedullary nails/bone marrow drills	•
Reprocessing of trays	•
Reprocessing of mesh trays	•
Reprocessing of measuring cylinders	•
Reprocessing of containers, e.g., soda-lime containers	•
Reprocessing of bedpans and urine bottles	•

APH 592
 Colour transfer belt
 for protocol printer APH 510 PRT110.



EAN: 8809166675907 / Material number: 11053500 / Old Material Number: 68817005D

Reprocessing of different utensils	•
Dimensions and weight	
External dimensions, net height in inches	0
External dimensions, net width in inches	0
External dimensions, net depth in inches	0
External dimensions, gross height in inches	2 1/4 (55)
External dimensions, gross width in inches	7 1/8 (180)
External dimensions, gross depth in inches	9 7/8 (250)
Net weight in lbs	1 (0.04)
Gross weight in lbs	1 (0.05)
Standards, test certificates, and symbols	
CE	•